



# UPLC, Sample Manager and QTOF System Performance Check - Proforma

---

## **NPC.PRO.MS002 Version 2.1**

---

Effective Date: April 2019

### **1. Purpose**

---

The purpose of this proforma is to document the Ultra Performance Liquid Chromatography (UPLC) Mass Spectrometer (MS) system instrument checks as outlined in the protocol NPC.SOP.MS002. The required sections of this proforma should be printed on the day of use, completed and then stored with all project specific documentation. This proforma should be used in combination with the SOP.

### **2. Proforma Approval**

---

---

Prepared by Dr Verena Horneffer-van der Sluis

---

Date

---

Approved by Dr Maria Gomez-Romero

---

Date

---

Authorised by Dr Matthew Lewis

---

Date

## UPLC and Q-TOF System Checklist

Please, indicate an action has been performed by ticking the appropriate box (☐)

Assay specifics			
Operator:		System No.:	
Project:		Date:	
LC system and Sample Manager – Visual inspection			
Non-essential material removed			<input type="checkbox"/>
LC system and Sample Manager have been checked and there are no leaks and/or corrosion			<input type="checkbox"/>
All solvent bottles are in date and topped up where necessary			<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>			
LC System and Sample Manager – Routine maintenance			
System (line A,B) primed for 4 minutes in methanol			<input type="checkbox"/>
Seal wash primed in parallel			<input type="checkbox"/>
Pressure stable ( $\Delta$ PSI <20 PSI)			Y/N
Dynamic leak test performed:			<input type="checkbox"/>
<b>Pump A</b>	Accumulator _____ nL/min	Primary _____ nL/min	Pass/Fail
<b>Pump B</b>	Accumulator _____ nL/min	Primary _____ nL/min	Pass/Fail
System (line A,B) primed for 4 minutes in fresh isopropanol			<input type="checkbox"/>
Weak and strong needle washes rinsed for 30 seconds			<input type="checkbox"/>
Syringe and sample loop washed			<input type="checkbox"/>
Sample Manager temperature checked and recorded _____ °C			<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>			

*Table continues*

Assay specifics			
Operator:		System No.:	
Project:		Date:	
MS system – Visual inspection			
System checked for leaks, vapour trap bottle checked and volume of LeuEnk solution at sufficient level			<input type="checkbox"/>
Mass spectrometer in operate mode and correct tune file loaded			<input type="checkbox"/>
Tune file _____			
MS tune page settings correspond to correct NPC standard settings for method			<input type="checkbox"/>
Method (RPOS/RNEG/HPOS/LPOS/LNEG): _____			
Vacuum read back pressure values			
Backing:	Collision:	TOF:	
Vacuum read back values acceptable?			Y/N
Comments: n/a <input type="checkbox"/>			
MS System – Check of LeuEnk solution			
Check if volume of LeuEnk solution is sufficient (>100 mL) and not expired? <i>If yes, go to next section. If not, prepare fresh as given in SOP NPC.SOP.MS001</i>			Y/N
Add 7.5 mL of LCMS grade water to the Waters LeuEnk bottle (3 mg) and mix until dissolved (400 ng/μL LeuEnk stock solution)			<input type="checkbox"/>
<input type="checkbox"/> n/a if stock already prepared. Prep date: _____			
Prepare 500 mL of 1:1 acetonitrile/H <sub>2</sub> O using a measuring cylinder (With/without 0.1% formic acid)			<input type="checkbox"/>
Add respective volume of 400 ng/μL LeuEnk stock solution using an automatic pipette			
Volume of 400 ng/ μL Leu Enk stock solution added:			<input type="checkbox"/>
Mix until the solution is homogenous and sonicate for 5 minutes			<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>			

Table continues

<b>MS System – Detector and Lockmass (<i>in assay specific polarity only</i>)</b>		
Fluidics line B purged with LeuEnk		<input type="checkbox"/>
<b>Detector Set-up (@ 5e<sup>5</sup> intensity)</b>		
	<i>Positive mode</i>	<i>Negative mode</i>
<b>Detector Voltage (V)</b>		
<b>Average Ion Area</b>		
Instrument uncalibrated		<input type="checkbox"/>
<b>Nominal mass</b> of LeuEnk model peak acceptable?		Y/N
<b>Veff corrected?</b>		Y/N
Stable LeuEnk signal intensity at 2e <sup>4</sup> in positive / negative mode		<input type="checkbox"/>
<b>Leu Enk Mass</b>		
	<i>Positive mode (556.2771 m/z)</i>	<i>Negative mode (554.2614 m/z)</i>
<b>Lockmass set up</b>		
	<i>Positive mode</i>	<i>Negative mode</i>
<b>Flow rate (μL/min)</b>		
<b>LockMass capillary volt (kV)</b>		
<b>DRE %</b>		
<b>Time (s)</b>		
<b>Comments:</b> n/a <input type="checkbox"/>		

**Table continues**

**MS system –Resolution check (*in assay specific polarity only*)**

Fluidics line A purged twice with SSTM (RPOS/RNEG/HPOS) or line B with LeuEnk (LPOS/LNEG)



Infuse at 15  $\mu\text{L}/\text{min}$  and ensure stable signal. Perform acquisition for 2.5 mins. Use ResCal to calculate the resolution

**Positive mode**

<b>Intensity</b>	$2e^4$	$2e^5$
<b>Flow (<math>\mu\text{L}/\text{min}</math>)</b>		
<b>Capillary voltage (kV)</b>		
<b>Res_1</b>		
<b>Res_2</b>		
<b>Res_3</b>		
<b>Res_4</b>		
<b>Res_5</b>		
<b>Res_average</b>		
Res_average value of $2e^4$ entered into Acquisition Settings		<input type="checkbox"/>

**Negative mode**

<b>Intensity</b>	$2e^4$	$2e^5$
<b>Flow (<math>\mu\text{L}/\text{min}</math>)</b>		
<b>Capillary voltage (kV)</b>		
<b>Res_1</b>		
<b>Res_2</b>		
<b>Res_3</b>		
<b>Res_4</b>		
<b>Res_5</b>		
<b>Res_average</b>		
Res_average value of $2e^4$ entered into Acquisition Settings		<input type="checkbox"/>

**Comments:** n/a ☐

**Table continues**

<b>MS system – Calibration and last checks (<i>in assay specific polarity only</i>)</b>		
Fluidics line C purged twice with sodium formate solution		<input type="checkbox"/>
<b>Calibration</b>		
	<i>Positive mode</i>	<i>Negative mode</i>
<b><i>RSM (ppm)</i></b>		
<b><i>95% Confidence band (ppm)</i></b>		
<b><i>Number of peaks acceptable</i></b>		
Purge lines A and/or C with wash solution twice and flush to remove residual sodium formate and SSTM		<input type="checkbox"/>
Back-up and save system settings		<input type="checkbox"/>
Check MS setting match assays specific values		<input type="checkbox"/>
Save tune file TOF___.ipr		<input type="checkbox"/>
Rest ADC-Results.csv		<input type="checkbox"/>
<b>Comments:</b> n/a <input type="checkbox"/>		