



RPC UPLC-QTOF Analysis of Lipids in Human Plasma and Serum - Proforma

NPC.PRO.MS003 Version 2.1

Effective Date: April 2019

1. Purpose

The purpose of this proforma is to document the Ultra Performance Liquid Chromatography (UPLC) Mass Spectrometer (MS) system reversed phase chromatographic (RPC) Lipid assay as outlined in the protocol NPC.SOP.MS003. This proforma should be used in combination with this SOP. The required sections of this proforma should be printed on the day of use, completed and then stored with all project specific documentation.

2. Proforma Approval

Prepared by Dr Verena Horneffer-van der Sluis

Date

Approved by Dr Maria Gomez-Romero

Date

Authorised by Dr Matthew Lewis

Date

Reagents

Assay specifics			
Analyst:			
Project:		Date:	
Chemical	Supplier, P/N	Batch/lot no.	Opened date
LPC(9:0)	Avanti, 855276P		
PC(11:0/11:0)	Avanti, 850330P		
FA(17:0)	Sigma, H3500		
PG(15:0/15:0)	Avanti, 840446P		
PE(15:0/15:0)	Avanti, 850704P		
PS(17:0/17:0)	Avanti, 840028P		
PA(17:0/17:0)	Avanti, 830856P		
Cer(d18:1/17:0)	Avanti, 860517P		
DG(19:0/19:0)	Sigma, 68633		
PC(23:0/23:0)	Avanti, 850372P		
TG(15:0/15:0/15:0)	Sigma, T4257		
TG(17:0/17:0/17:0)	Sigma, T2151		
LCMS grade acetic acid			
LCMS grade ammonium acetate			
LCMS grade acetonitrile			
LCMS grade isopropanol			
LCMS grade water			
Phosphoric acid (85-90%)			

PART A - Lipid Standard Mixture Preparation

Please, indicate an action has been performed by ticking the appropriate box (☐)

Assay specifics					
Analyst:					
Project:				Date:	
Lipid stock solutions					
Weigh each lipid standard individually into the corresponding volumetric flask outlined below and record the weight and volume (quantities can be scaled up or down depending on requirements)					<input type="checkbox"/>
Make each volumetric flask up to volume with isopropanol					<input type="checkbox"/>
Lipid	Weight (mg)	Volume of isopropanol (mL)	Target stock conc. (mg/mL)	Weight (mg)	Volume (mL)
LPC(9:0)	5	10	0.5		
PC(11:0/11:0)	5	50	0.1		
FA(17:0)	10	100	0.1		
PG(15:0/15:0)	5	500	0.01		
PE(15:0/15:0)	5	500	0.01		
PS(17:0/17:0)	20	200	0.1		
PA(17:0/17:0)	5	50	0.1		
Cer(d18:1/17:0)	5	50	0.1		
DG(19:0/19:0)	40	200	0.2		
PC(23:0/23:0)	5	50	0.1		
TG(15:0/15:0/15:0)	10	100	0.1		
TG(17:0/17:0/17:0)	10	100	0.1		
Sonicate stock solutions until sample dissolution observed					<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>					

Assay specifics

Analyst:				
Project:			Date:	
Lipid standard mixture				
Transfer each lipid standard stock solution prepared above to a 2 L volumetric flask according to the quantities outlined in the table below using appropriate automatic pipettes or measuring cylinders				
Lipid	Vol. of stock in mix (mL)	Target conc. (µg/mL) in 2 L of isopropanol	Final conc. (µg/mL) in 2 L	
LPC(9:0)	1	0.25		<input type="checkbox"/>
PC(11:0/11:0)	5	0.25		<input type="checkbox"/>
FA(17:0)	50	2.50		<input type="checkbox"/>
PG(15:0/15:0)	200	1.00		<input type="checkbox"/>
PE(15:0/15:0)	50	0.25		<input type="checkbox"/>
PS(17:0/17:0)	120	6.00		<input type="checkbox"/>
PA(17:0/17:0)	20	1.00		<input type="checkbox"/>
Cer(d18:1/17:0)	1	0.05		<input type="checkbox"/>
DG(19:0/19:0)	120	12.00		<input type="checkbox"/>
PC(23:0/23:0)	5	0.25		<input type="checkbox"/>
TG(15:0/15:0/15:0)	50	2.50		<input type="checkbox"/>
TG(17:0/17:0/17:0)	50	2.50		<input type="checkbox"/>
Make the 2 L volumetric flask up to volume with isopropanol				<input type="checkbox"/>
Mix the volumetric flask until the contents is visually homogenous				<input type="checkbox"/>
Aliquot the Lipid standard mixture and store at -20 °C until required				<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>				

PART B - Long Term Reference (LTR), Study Reference (SR) and SR Dilution Series Preparation - *On the Day of MS-SR preparation*

Please, indicate an action has been performed by ticking the appropriate box (☐)

Assay specifics	
Analyst:	
Project:	Date:
LTR	
Remove per sample batch one 11.5 mL stock plasma/serum LTR from storage at -80 °C and allow to defrost for 2 hours at 2-8 °C (1 tube per 2 sample batches)	<input type="checkbox"/>
Aliquot 450 µL of plasma/serum LTR in 25 separate Eppendorf tubes	<input type="checkbox"/>
Store at -80 °C until required	<input type="checkbox"/>
SR	
Aliquot 450 µL of undiluted pooled MS-SR per sample plate in _____ separate Eppendorf tubes.	<input type="checkbox"/>
Aliquot 100 µL of SR into at least 7 separate Eppendorf tubes (1 set per sample batch, ≤ 1000 samples) for instrument conditioning, start/end SR, and DIDA	<input type="checkbox"/>
Store at -80 °C until required	<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>	

Table continues

Assay specifics

Analyst:

Project:

Date:

No. of samples: _____ No. of batches: _____ (Sample batch is ≤1000 samples)

Total number of dilution series sets required (sample batch x (4 + 2 backup)) : _____

The table below indicates the volumes required for a single sample set (6 dilution series). This will be followed for studies consisting of 1000 samples or less. Please cross through and include an appropriate table as an attachment to this document if the number of samples in the study is more than 1000 samples

For ≤1000 samples

Dilution point	Percentage of SR (%)	Vol. of SR (μL)	Vol. of LCMS grade water (μL)	Total vol. (μL)	Vol. in aliquot (μL)	
1	100	450	0	450	70	<input type="checkbox"/>
2	80	240	60	300	45	<input type="checkbox"/>
3	60	150	100	250	35	<input type="checkbox"/>
4	40	100	150	250	35	<input type="checkbox"/>
5	20	100	400	500	45	<input type="checkbox"/>
6	10	50	450	500	70	<input type="checkbox"/>
7	1	10	990	1000	70	<input type="checkbox"/>
Blank	0	0	500	500	70	<input type="checkbox"/>

Aliquot each dilution point into 6 separate Eppendorf tubes according to the volumes detailed above (Vol. in aliquot) and store at -80 °C until required.

☐

Comments: n/a ☐

PART C - SR Dilution Series, SR and Blanks - *On the Day of Analysis*

Please, indicate an action has been performed by ticking the appropriate box (☐)

Assay specifics				
Analyst:				
Project:			Date:	
Remove 1 set of SR dilution series aliquots per polarity from storage at -80 °C, appropriate instrument conditioning or DIDA SR, and blanks from storage at -80°C and allow to defrost for 2 hours at 2-8 °C				<input type="checkbox"/>
Remove Lipid IPA standard mix from storage -20 °C and allow to defrost for 2 hours at 2-8 °C.				<input type="checkbox"/>
Dilution Series				
Complete all pre-project system performance checks as given NPC.SOP.MS002				<input type="checkbox"/>
Add the corresponding volume of the lipid standard mixture to the SR dilution aliquot as per the table below using pipettes, then vortex mix				<input type="checkbox"/>
Dilution point	Percentage of SR (%)	Vol. in aliquot (µL)	Vol. of Lipid mix (µL)	
1	100	70	280	<input type="checkbox"/>
2	80	45	180	<input type="checkbox"/>
3	60	35	140	<input type="checkbox"/>
4	40	35	140	<input type="checkbox"/>
5	20	45	180	<input type="checkbox"/>
6	10	70	280	<input type="checkbox"/>
7	1	70	280	<input type="checkbox"/>
Blank	0	70	280	<input type="checkbox"/>
Mix for 2 hours on a plate mixer at 1400 rpm at 4 °C				<input type="checkbox"/>
Centrifuge at 3486 g for 10 minutes at 4°C				<input type="checkbox"/>
Aliquot supernatant into labelled 300 µL vials				<input type="checkbox"/>
Place vials in correct positions in autosampler				<input type="checkbox"/>

Table continues

Assay specifics	
Analyst:	
Project:	Date:
100 µL SR aliquot(s) and Lipid IPA were allowed to defrost for 2 hours at 2-8 °C. Number of Eppendorf tube defrosted: _____	
<input type="checkbox"/>	
Instrument Conditioning or DIDA (100% sample)	
Add 400 µL of Lipid IPA standard mixture in each tube	<input type="checkbox"/>
Mix for 2 hours on a plate mixer at 1400 rpm at 2-8 °C	<input type="checkbox"/>
Centrifuge at 3486 g for 10 minutes at 4 °C	<input type="checkbox"/>
Pipette supernatant into labelled 300 µL vials and seal with foil lined caps	<input type="checkbox"/>
Start/End SR (50% sample)	
Add 100 µL of LCMS grade water to each tube	<input type="checkbox"/>
Mix for 5 min at 1400 rpm at 2-8 °C	<input type="checkbox"/>
Add 100 µL of Lipid IPA standard mixture to each tube	<input type="checkbox"/>
Mix for 2 hours on a plate mixer at 1400 rpm at 2-8 °C	<input type="checkbox"/>
Centrifuge at 3486 g for 10 minutes at 4 °C	<input type="checkbox"/>
Pipette supernatant into labelled 300 µL vials and seal with foil lined caps	<input type="checkbox"/>
Blanks	
Pipette 100 µL of LCMS grade water in an Eppendorf	<input type="checkbox"/>
Add 400 µL of Lipid IPA standard mixture	<input type="checkbox"/>
Mix for 2 hours on a plate mixer at 1400 rpm at 2-8 °C	<input type="checkbox"/>
Centrifuge at 3486 g for 10 minutes at 4 °C	<input type="checkbox"/>
Pipette into labelled 300 µL vials and seal with foil lined caps	<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>	

PART D - Sample Preparation

Please, indicate an action has been performed by ticking the appropriate box (☐)

Assay specifics		
Analyst:		Date:
Project:		Plate Identity:
1	Thaw sample plate, LTR, SR and Lipid IPA standard mixture at 2-8 °C for of 2 hours	<input type="checkbox"/>
2	Centrifuge lipid plate at 3486 g for 1 minute at 2-8 °C	<input type="checkbox"/>
3	Remove carefully seal cap mat and place side up on a clean surface	<input type="checkbox"/>
4	Add 50 µL of plasma/serum LTR into column 11 of the plate	<input type="checkbox"/>
5	Add 50 µL of plasma/serum SR into column 12 of the plate	<input type="checkbox"/>
6	Add 50 µL of LCMS grade water into each well	<input type="checkbox"/>
7	Cover the plate with the seal cap mat	<input type="checkbox"/>
8	Mix for 5 min on a plate mixer at 1400 rpm at 2-8 °C	<input type="checkbox"/>
9	Centrifuge lipid plate at 3486 g for 1 minute at 2-8 °C	<input type="checkbox"/>
10	Remove carefully seal cap mat and discard	<input type="checkbox"/>
11	Add 400 µL of the Lipid IPA standard mixture to each well	<input type="checkbox"/>
12	Heat seal the sample plate with foil seal	<input type="checkbox"/>
13	Mix for 2 hours on a plate mixer at 1400 rpm at 2-8 °C	<input type="checkbox"/>
14	Centrifuge sample plate at 3486 g for 10 minutes at 4 °C	<input type="checkbox"/>
15	Label three 96 well microplates, one for positive, one for negative and one for backup; with LPOS, LNEG, BACKUP, project name, date and barcode label	<input type="checkbox"/>

Table continues

Continuation -PART D - Sample Preparation

16	Carefully remove heat seal foil from the sample plate (<i>without disturbing the pelleted material</i>)	<input type="checkbox"/>
17	Transfer 100 µL of each sample to the corresponding well in the analytical plates (LPOS, LNEG and BACKUP)	<input type="checkbox"/>
18	Heat seal all analytical plates with heat seal foil	<input type="checkbox"/>
19	Log the plates in LIMS	<input type="checkbox"/>
20	Centrifuge LPOS and LNEG plates at 3486 g for 5 minutes at 4 °C	<input type="checkbox"/>
21	Store BACKUP plate at -80 °C	<input type="checkbox"/>
22	Place LPOS and LNEG plates in relevant autosamplers	<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>		

PART E – Preparation of Mobile Phases and Wash Solutions

Please, indicate an action has been performed by ticking the appropriate box (☐)

Assay specifics	
Analyst:	
Project:	Date:
Preparation of lipid mobile phase	
N.B. Volumes can be scaled up or down depending on requirement	
500 mM phosphoric acid solution	
Transfer 10 mL of LCMS grade water into a Duran bottle Actual Volume = _____	<input type="checkbox"/>
Add 340 µL of 85 – 90% phosphoric acid using a pipette to the Duran bottle from above Actual Volume = _____	<input type="checkbox"/>
Mix until the content is homogenous, seal the bottle and assign an expiry of 1 month	<input type="checkbox"/>
Mobile phase A: preparation of 2 litres (water:isopropanol:acetonitrile 2:1:1, 5 mM ammonium acetate, 0.05% acetic acid, 20 µM phosphoric acid)	
Weigh 0.7708 g ± 0.0077 g of ammonium acetate into a weigh boat Weight = _____	<input type="checkbox"/>
Transfer the ammonium acetate into a 2 L Duran bottle	<input type="checkbox"/>
Using a measuring cylinder, transfer 1 L of LCMS grade water to the Duran bottle	<input type="checkbox"/>
Sonicate for 5 minutes, or until the ammonium acetate has fully dissolved	<input type="checkbox"/>
Add 1 mL of acetic acid to the Duran bottle using a Hamilton syringe	<input type="checkbox"/>
Add 80 µL of the 500 mM phosphoric acid solution to the Duran bottle and mix until homogenous	<input type="checkbox"/>
Measure the pH of the mobile phase ensuring the pH is 4.4 ± 0.1 pH = _____	<input type="checkbox"/>
Using a volumetric flask, slowly transfer 500 mL of isopropanol to the Duran bottle (with mixing)	<input type="checkbox"/>
Using a volumetric flask, slowly transfer 500 mL of acetonitrile to the Duran bottle (with mixing)	<input type="checkbox"/>
Mix until the solution is completely homogenous	<input type="checkbox"/>
Sonicate for 10 minutes	<input type="checkbox"/>

Table continues

Assay specifics	
Analyst:	
Project:	Date:
Mobile phase B: preparation of 2 litres (<i>isopropanol:acetonitrile 1:1, 5 mM ammonium acetate, 0.05% acetic acid</i>)	
Weigh 0.7708 g \pm 0.0077 g of ammonium acetate into a weigh boat Weight = _____	<input type="checkbox"/>
Transfer the ammonium acetate into a 2 L Duran bottle	<input type="checkbox"/>
Using a measuring cylinder, transfer 1 L of LCMS grade isopropanol to the Duran bottle	<input type="checkbox"/>
Sonicate for 45 minutes, or until the ammonium acetate has fully dissolved	<input type="checkbox"/>
Add 1 mL of acetic acid to the Duran bottle using a Hamilton syringe	<input type="checkbox"/>
Using a volumetric flask, slowly transfer 1 L of acetonitrile to the Duran bottle (with mixing)	<input type="checkbox"/>
Mix until the solution is completely homogenous	<input type="checkbox"/>
Sonicate for 10 minutes	<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>	

Table continues

Assay specifics

Analyst:

Project:

Date:

Preparation of lipid wash solutions

N.B. Volumes can be scaled up or down depending on requirement

Seal wash (*isopropanol:water 1:9*)

Transfer 100 mL of isopropanol into a Duran bottle

☐

Add 900 mL of LCMS grade water into the Duran bottle

☐

Mix until the content is homogenous

☐

Sonicate for 5 minutes, seal the bottle and assign an expiry of 1 month

☐

Weak needle wash (*water:isopropanol 1:4*)

Transfer 200 mL of LCMS grade water into a Duran bottle

☐

Add 800 mL of LCMS grade isopropanol into the Duran bottle

☐

Mix until the content is homogenous

☐

Sonicate for 5 minutes, seal the bottle and assign an expiry of 1 month

☐

Strong needle wash (*isopropanol*)

Use isopropanol as supplied and assign an expiry of 3 months

☐

Comments: n/a ☐

PART F - Acquisition

Please, indicate an action has been performed by ticking the appropriate box (☐)

Assay specifics	
Analyst: _____	
Project: _____	Date: _____
Instrument number: _____	
Column: Acquity UPLC BEH C8 1.7µm, 2.1x100 mm; P/N:186002878 LOT: _____ Serial Number: _____	
Ionisation mode required: _____	
Instrument check performed (see separate Proforma sheet NPC.PRO.MS002)	<input type="checkbox"/>
All solvent lines match the assay specific buffers and solutions	<input type="checkbox"/>
Check the following against NPC.SOP.MS003:	
Tune file used: _____	
Acquisition method used: _____	
Inlet method used: _____	
Sequence loaded into Masslynx	<input type="checkbox"/>
Sample plate loaded into tray position	<input type="checkbox"/>
Comments: n/a <input type="checkbox"/>	